

# PERFORMANCE

## 2 & 3 SECTION



Our robust entry level Performance range of treatment couches are available in either 2 or 3 sections with electric height operation, making them ideal for a wide range of healthcare environments and treatments procedures.

**2 Section: TMP2E** - Electric height elevation    **3 Section: TMP3E** - Electric height elevation

### STANDARD SPECIFICATION

- Safe working load and lifting capacity of 170kg (26 stone)
- Excellent height range from 50cm (20") to 95cm (37")
- Hoist clearance along entire length of couch
- 125mm diameter individually locking castors
- 50mm high density flame-retardant foam
- Antimicrobial upholstery
- Fully compliant and CE marked
- Hand switch fitted as standard
- Gas spring assisted back rest from -30° to +90°



### OPTIONAL ACCESSORIES

- Single foot switch operation
- Paper towel holder
- Phlebotomy arm
- Breathing hole

Safe working load and lifting capacity of 170kg (26 stone)

### DIMENSIONS

- Length: 190cm
- Width: 66cm
- Safe Load: 170kg
- Height range: 50cm-95cm




### WARRANTY

- 1 year warranty on moving components
- 3 year warranty on frame

### STANDARDS

- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- IP rating - IPx4
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

#### Performance 2 Section:



#### Performance 3 Section:

