

# ULTRASOUND MEDICAL COUCH



An electric height and motorised three section ultrasound couch featuring a Trendelenburg tilt to allow for patient scanning. The couch features 55mm retractable castors with locking function as well as a 90° drop for patient access.

## TMUS - Ultrasound Medical Couch

### STANDARD SPECIFICATION

- Electric height, motorised backrest and trendelenburg tilt adjustment
- Leg section adjustable to -90°, seat tilt from 0° to +25° & backrest adjustment from -15° to +90°
- Stay level swing back arm rests for easy patient access (arms also removable)
- Leg length adjustment +15cm (6") on divided leg model
- Couch goes horizontal for patient recovery in emergency situations
- Purpose made precision bearings on all pivot joints
- Footrest adjustment by gas strut mechanism
- Adjustable foot for uneven floors
- Retractable castor wheels with locking function

### OPTIONAL ACCESSORIES

- Breathing hole and plug
- Paper roll holder
- Double thickness padded upholstery
- Matching operators chair or stool

### STANDARDS


- IP rating - IPx4
- EEC classification - MDD93/42/EEC Annex IX Class I
- Electrical standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl flame retardant to CRIB5 (BS7176 inc. foam)
- Foam to BS 5852: Part 2: 1982, Ignition Source 5 (CRIB5)

### WARRANTY


- 3 year warranty on moving components
- Lifetime warranty on frame
- Fully compliant and CE marked

### DIMENSIONS

- Length: 174cm, Width: 63cm
- Height range: 45cm - 106cm
- Safe Load: 250kg



Height Range 45cm - 106cm



Width 63cm



Safe working  
load and lifting  
capacity of 250kg  
(39 stone)



Telephone: +44 (0) 1484 865586, Email: [info@titanmed.co.uk](mailto:info@titanmed.co.uk), Web: [www.titanmed.co.uk](http://www.titanmed.co.uk)  
TitanMed - part of The Steri-Pod Group, Unit 8, Sovereign Business Park, Huddersfield HD8 8FW  
TitanMed Limited Registered in England 10120881, VAT Registration Number 272566288